## Accessories

# MEMPROTEC



#### Syrup filter press 17x17

Several sizes available up to 10 filtration plates with hydraulic jack and 12 ton air.

## Syrup filter press 20x20

Several sizes available from 7 to 20 filtration plates with hydraulic jack and 20 ton air.

#### Accessories included with the presses:

- Drain valve
- Industrial wheels on bearing
- Tilt failure for better emptying
- Drain plug on failure
- Stainless steel jack rod
- Knife to reduce number of plates
- 2 sets of filters for each syrup press
- Pump mounting bracket

## Complete filter kit glass bead media



| Plastic filter 45 gpm         | 399954  |
|-------------------------------|---------|
| Plastic filter 75 gpm         | 399955  |
| Stainless steel filter 75 gpm | 399957  |
| Glass media 60 lbs            | 3807660 |

The glass media filter allows maple sap to be filtered before it enters in the Reverse Osmosis System. This ensures well-filtered sap and helps reduce the amount of Purtrex filter cartridges needed for filtration. It is still important to always use the 5 micron 20" filter cartridges on the Reverse Osmosis System to ensure good filtration.

The filter is available in plastic or stainless steel. The pumps are also available in plastic or stainless steel.

Memprotec glass bead media blend provides superior filtration by retaining particles 3 microns and larger. The filter is easy to backwash in minutes with a single, simple, easy-to-use control valve. On the stainless-steel model, there is a 5" porthole at the top to fill the glass bead media filter as well as a second porthole at the bottom to empty it. It also has a drain valve with strainer for drainage.



### **Cleaners and accessories**

Monitoring and maintenance of your membranes every day of the sugar season is essential to maintain optimal concentration performance throughout the season and to preserve your membranes in the long run. **MEMPROTEC** offers you a complete range of high performance maintenance products for your RO.

### Membrane 8"

| Membrane MemPro HC+8040     | 299701  |
|-----------------------------|---------|
| Membrane MemPro HC-70       | 299686  |
| Membrane MemPro HP-70       | 299700  |
| Membrane MemPro 70-400      | 3999915 |
| Membrane MemPro 8040-LFT-HE | 399991  |
| Membrane MemPro SF-70       | 2996996 |

#### Membrane 4"

 Membrane HP-40
 299687

 Membrane MemPro 4040-LFT-HE
 399994

 Membrane MemPro SF-40
 2996997

| Filter cartridge     |                     |        |       |
|----------------------|---------------------|--------|-------|
| K                    | 9 -7/8" Purtrex 5 m | icrons | 39905 |
| L                    | 20" Purtrex 5 micro | ons    | 38002 |
| 4½ x 20" Soap filter |                     |        |       |
| 0                    | Single use          | 255496 |       |
| J                    | <b>Reusable SS</b>  | 380449 |       |

#### 7" x 30" Soap prefilter

00

380775





| Cleaners for membranes   |               |         |      |      |         |
|--|---------------|---------|------|------|---------|
| Bio-Clean B cleaner  |               |         |      |      |         |
| Membrane cleaner for Mempro 8040 LFT-HE, Filmtec,<br>Mempro NF-70, H20, Hydranautics |               |         |      |      |         |
| A  | 4 L           | 380015  | B    | 20 L | 380016  |
| Bio  | Bio-Protec PG |         |      |      |         |
| Long term storing solution.  |               |         |      |      |         |
| C  | 4 L           | 37985PP | D    | 20 L | 37986PP |
| Rio  | -Acid         |         |      |      |         |
| Bio-Acid   |               |         |      |      |         |
| Cleaner for reverse osmosis or nanofiltration membrane.                              |               |         |      |      |         |
| E  | 2 KG          | 37995   | LIQU | ID   |         |
| F  | 4 KG          | 379955  | 20 L | 380  | 515     |

205 L

380514

#### **Bio-Preserve**

18 KG

G

379951

Short term preservation.

H 16 oz 379905

#### **OXISAN Cleaner**

**4 L** 38006 **20 L** 380061

### **Cleaners for pans**





### **Cleaners for tubing**

### Isopropyl Alcohol 70%

**20 L** 380340 **205 L** 380342



## Accessories



| Measuring tools   | Inline refractometer                             | Automatic draw-off valve  |
|---|--|---|
| Digital Clamp Meter 8344707   | <b>0-33</b> 3603                                 | 506050  |
| Digital refractometer 694076  | <b>0-80</b> 3604 (for syrup)                     |   |
| 11 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1  | 3603 3604  |   |
| BARINOX Stainless Steal Drums   | Strainer   | Valve panel   |
| <b>25 gl</b> 3805100<br><b>34 gl</b> 3805101<br><b>40 gl</b> 3805104  | 1" 398505<br>1½" 39957<br>2" 398504<br>3" 398507 | 1/2"       501525         3/4"       501526         1"       501527         1-1/4"       501530 |
| Barrel cap 3805102  |  | Mombuser charge Voccol  |
| 20" Filter Housing<br>Big blue 399931   | Electrical components                            | Membrane storage Vessel   |
| <b>10</b> " 38019<br><b>20</b> " 38009  |  | <b>4''</b> 899254<br><b>8''</b> 899255  |
|   | GFCI 240 volts 20 amp (MemPro 350) 2021605       |   |
|   | 120V Dual Wall GFCI outlet200406                 |   |
| Membrane tools  | Memproflush regulator                            | SS Filter   |
| Membrane puller 8"       693862         Membrane puller 8 ½"       693863         Membrane pusher 8 ½"       693863 | 399575M  | SS Filter with 7 cartridges of 20"<br>501910  |